

## F423 --- Enter & Store a SimpleVital (PSI 8)

F414_US8307_ F423 - Enter & Store a Simple Vitals			
<b>Test Title:</b> Enter & Store a Simple Vital <b>Description:</b> This covers the functionality of Enter & Store a Simple Vital.			
<b>Pre-Condition(s):</b> User must have access to UI-components. User must have CPRS account with eHMP write up role.			
<b>Post-Condition(s):</b> N/A			
Step Title	Description	Expected Results	Data
Access the Observations Tab	After logging into eHMP and selecting/confirming a patient, select the Observations tab from the Patient Care Information component.	A form is displayed with the heading, "All Observations".  (Note : "Delete" and "view" on the button of the footer is enabled only when one items from the list are selected.)	
Create New Observation	User selects the "Create New Observation" button in All Observations tray.	A "Create New Observation" tray header opens with a list of options for selection. Make sure all the fields are according to mock up <a href="http://10.1.1.150/ui-components/views/psi8/store-vital-tray.html">http://10.1.1.150/ui-components/views/psi8/store-vital-tray.html</a> Add new observation (heading) Vitals Immunization Meds, Non-va Problems Vitals	

Access Vitals Form	From the List of Items user selects "Vitals".	The "Enter Vitals" forms appear. Make sure all the fields are according to mock up.	
Pass button	User selects "Pass" button	All the fields on the form are disabled and the button changes to "Deselect All" button	
Deselect All	User selects "Deselect All" button	All the fields on the form are "Enabled" and the button changes to "Pass" button	
Expand All Icon	User selects one of the "Expand all" icons on the form.	The appropriate section is expanded to show all of the fields and the button is changed to "Collapse All"	
Collapse All Button	User selects the Collapse All button.	The grey fields within each section are hidden from the view and the button is changed to "Expand All"	
Footer	At the footer, user can select one of the buttons (Delete, Close and Add).	All buttons appear within the footer of the "Vitals" form.	
Delete	User selects the "Delete" button from the footer.	Alert form with header and text disabled with a option to "cancel" button and "continue" button	
Cancel Alert	Repeat step 7 and then User selects the "cancel" button from the footer.	Focus should go back to the Cancel button on Entered in Error Form	

Add button	User repeats steps 1-3 and enters data into the form. User selects the Add button.	A growl alert with a confirmation message appears.  (Note: Make sure all the fields are exist according to wireframe and the focus goes to "Vitals Added" growl alert.)	
"X" button	User repeats steps 1-3 and selects the "X" button on top of the form are selects	The form is closed and focus goes back to tray.	
Mandatory fields	Within the "Vitals" tray make sure Positive and negative steps for mandatory fields are completed	If user does not enter a mandatory field and selects "add" or "save" the fields should be highlight in red with a error message in the bottom	
. Options Icon	Select the Options icon on top right corner of the form.	The list of options for the note appear in a dropdown (Preview, Print, Delete)	