

F432 --- Enter & Store a SimpleLab orders (PSI8)

F432 --- Enter & Store a Simple Lab orders			
Test Title: Enter & Store a Simple Lab orders Description: This covers the functionality of Enter & Store a Simple Lab orders			
Pre-Condition(s): User must have access to UI-components. User must have CPRS account with eHMP write up role.			
Post-Condition(s): N/A			
Step Title	Description	Expected Results	Data
1. Access the Record and Enter & Store a Simple Lab orders	After logging into eHMP and selecting/confirming a patient, select the Observations tab from the Patient Care Information component.	A form is displayed with the heading, "All Observations". (Note : "delete" and "view" on the button of the footer is enabled only when one items from the list are selected)	
2. New Entry	User selects the "New Entry" button in All Observations tray.	A "New Entry" tray header opens with a list of options for selection. Make sure all the fields are according to mock up Add new observation (heading) Lab orders Immunization Meds,Non-va Problems Lab orders	
3. Selection "Lab orders"	From the List of Items user selects "Lab orders".	A "Enter Lab orders" tray header appears. Make sure all the fields are according to mock up.	
4. "Back to Order Summary" button	User select "Back to Order Summary" button from footer of "Add a New Observation" form.	Focus goes back to "All Observations" Tray header.	

5. Mandatory fields	Within the "Lab orders" tray make sure Positive and negative steps for mandatory fields are completed	If user does not enter a mandatory field and selects "Accept" the fields should be highlight in red with a error message.	
6. "Cancel" button footer	User selects "Cancel" button from footer.	Conformation cancellation form will appear, User can select "cancel" or "continue"	
7. Cancel Alert	Repeat step 6 , User selects "cancel"	Focus goes back to vitals form footer	
8. Continue Alert	Repeat step 6, User selects "continue"	A growl alert will appear with a confirmation message.	
9. "Accept" button	User selects "Accept" drop down button.	There will be two dropdown options for the user to select from "Accept and add other" and "Accept and "sign"	
10. "Accept" button	Repeat step 9, User selects "Accept and add other" option from the drop down	The form will be add and a new Lab order form will be open	

11. "Accept" button	Repeat step 9, "Accept and "sign" option from the drop down	The form will be saved and sign task Form will be open.	
12. Close "X" button	User selects the Close or "X" button on top of the form.	A growl alert with "Close and Save" confirmation message will appear. Make sure all the fields are exist according to wireframe and then the user is returned to the "All Observations" Tray header.	

Observations Summary Lab orders Details View			
13. Observations Summary Page	Observations Summary Page will contain all of the Observations content sorted by Status (Default view). Meaning all drafts (Lab orders, Problem, Lab orders) within the Observations Summary are combined and selectable.	Make sure all the fields are according to mock up	
14. Arrange by: status	User clicks on "Arrange by :status" drop-down button	User can sort by 4 option "Location" "Provider" "status" and "title"	
15. Most Recent	User clicks on "Most Recent " drop-down button	Most Recent will change to "least recent " Note : user can sort by select "Most Recent " or "least recent	
16. "Lab orders" entry	Instead of selecting "New Entry " button on "All observations" tray user selects "Lab orders" entry.	Lab orders Entry List will open.	
17. "delete" buttons	Selection of a certain type of record will provide the following buttons in the Summary View below: "Delete", "edit" and "Entered in Error" when user selects delete	If a user opens a form and they never enter any data, but selected Delete. The form will close and return the user back to the button that initiated the form. If the user entered information and it is in draft mode, but selected Delete. A form(Growl	

		Alert" will appear asking "Are you sure you want to delete"? Yes or No	
18. "X" button	The user selects the "X" button.	Form closes and focus goes back to back to the button that initiated the form.	
19. "Edit" button	On the action bar user select "Edit" button	A "Lab orders" tray header appears. User can now make the changes.	